Micro I Diagnostic Spirometer

Operating Manual – English

Operating Manual – CareFusion USA

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Contents

I.	Introduction – English	9
	Package contents	10
	Contraindications, Warnings and Cautions	11
	Intended Use	12
	Environment	12
	Getting Started	12
	Configuration	13
	Micro I PC Software	14
	Operation	18
	Main Menu Overview	18
	Quick Exam	19
	Exam with Predicted values	21
	Post BD Exam	24
	NLHEP Mode	24
	NLHEP Quality Checks	25
	NLHEP QC Grades	26
	NLHEP Interpretation	26
	Switching Off	27
	Maintenance	27
	Battery Management	27
	Battery Replacement	28
	Calibration Check	30
	Cleaning the Transducer	31
	Servicing	31
	Trouble Shooting Information	32
	Electromagnetic Compatibility (EMC)	32
	Safety Designation per IEC60601-1	33
	Symbols	34
	Specification of the Micro I	36
	Accessories	37
	Customer Contact Information	212

Contents

II. Introduction – CareFusion USA	38
Package contents	39
Contraindications, Warnings and Cautions	40
Intended Use.....	41
Environment.....	41
Getting Started.....	41
Configuration	42
Micro I PC Software.....	43
Operation	47
Main Menu Overview	47
Quick Exam.....	48
Exam with Predicted values.....	50
Post BD Exam.....	53
NLHEP Mode.....	53
NLHEP Quality Checks	54
NLHEP QC Grades	55
NLHEP Interpretation.....	55
Switching Off.....	56
Maintenance.....	56
Battery Management	56
Battery Replacement.....	57
Calibration Check.....	59
Cleaning the Transducer	60
Servicing	60
Trouble Shooting Information	61
Electromagnetic Compatibility (EMC)	61
Safety Designation per IEC60601-1.....	62
Symbols	63
Specification of the Micro I	65
Accessories	66
Customer Contact Information.....	212

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Operating Manual – English

I. Introduction

The Micro I is a compact, rechargeable battery operated and fully portable diagnostic spirometer. It is accurate to the requirements of the ATS/ERS TASK FORCE: STANDARDISATION OF LUNG FUNCTION TESTING 2005. Its ergonomic and user-friendly design allows diagnostic spirometry measurements, including predicted values and automatic interpretation, to be made quickly and simply.

The Micro I features include:

- Measures FEV1, FVC(FEV6), FEV1/FVC, FEV1/FEV6, PEF, FEF25, FEF75, and FEF 25-75
- Fully configurable using software supplied
- ECCS, Asian or NHANES III predicted values
- Variation from norm as percentage or Z-score
- NLHEP, NICE, or the ATS/ERS interpretation
- Spirometry manoeuvre quality checks
- Post bronchodilator comparison
- NLHEP compliant mode

The spirometer uses the CareFusion Digital Volume Transducer, an extremely stable form of volume transducer, which measures expired air directly at B.T.P.S (Body Temperature and Pressure with Saturated water vapour) thus avoiding the inaccuracies of temperature corrections. This transducer is insensitive to the effects of condensation and temperature and avoids the need for individual calibration prior to performing a test.

CareFusion can supply spirometers to fulfil all your diagnostic and monitoring spirometry needs. The Micro II spirometer is the sister spirometer to the Micro I specifically designed for monitoring of asthma and COPD either at home or in the work place.

Package Contents

The Micro I is packaged in a convenient carrying case and comes complete with the following items (Fig.1):

1. Micro I spirometer
2. CareFusion Digital Volume Transducer together with disposable cardboard mouthpieces and instruction manual.
3. Universal power supply (36-PSU1017 5VDC 1.2A)
4. Cardboard mouthpiece adapter
5. USB/Charging cable

Fig 1

Contraindications, Warnings and Cautions

Contraindications: As performing spirometry might be physically demanding, it is contraindicated in patients with recent myocardial infarction. Also extensive exhalation might lead to syncope.

The following terms are used as follows in this manual

Caution: Possibility of injury or serious damage

Warning: conditions or practices that could result in personal injury.

Please Note: Important information for avoiding damage to the instrument or facilitating operation of the instrument.

	CAUTION: Read the manual before use
	WARNING: The instrument is not suitable for use in the presence of explosive or flammable gases, flammable anaesthetic mixtures or in oxygen rich environments.
	CAUTION: Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement.
	CAUTION: Pulmonary filters are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may also increase air resistance and lead to an incorrect measurement.
	CAUTION: Do not allow the patient to handle the spirometer when connected to either the power supply for charging or to a PC when configuring the unit.
	PLEASE NOTE: The product you have purchased should not be disposed of as unsorted waste. Please utilise your local  WEEE collection facilities for the disposal of this product.

Intended Use

The Micro I is used by clinicians and healthcare professionals for the diagnosis of asthma and COPD in adults and children over the age of 3 years.

Environment

Please observe the following precautions:

- Avoid exposing the Micro I to direct sunlight.
- Avoid operating the spirometer in dusty conditions or near to heating appliances or radiators.
- Do not keep the spirometer in a damp place or expose it to extreme temperatures.
- Do not direct the transducer holder towards a strong light source whilst operating the spirometer.

Getting Started

It is recommended that Micro I spirometer is fully charged before use. The power supply is provided with separate UK, USA and European plugs. Connect the required plug to the power supply and plug in to a mains outlet. Connect the Micro I to the power supply using the USB/charging cable and the charging symbol will appear on the screen.

Fully charge for a minimum of 5 hours when used for the first time.

Remove the protective film from the display screen before use.

Configuration

When the unit is first turned on it may be configured for your region. This will set the language, height and weight units, date format and predicted values set appropriate for your region. It will also configure the indices to be displayed and whether percentage predicted or Z-score is to be used to show variation from the norm. However all of these settings may be customised using the PC software supplied.

Turn the unit on by pressing the on/off button located at the top of the device and the following will be displayed:

Use the up and down arrow keys to highlight the required country and then press 'Enter' (↵). This procedure is only required when the unit is first switched on and the selection will be stored for future use.

It is recommended that the PC software be used to adjust these settings, if required. Connect Micro I to the PC using the USB cable supplied. Run the software, and turn on the Micro I. The PC software will detect that the Micro I is connected. Follow the on-screen instructions to configure your unit.

Micro I PC Software

This software allows you to customise the settings on your Micro I spirometer. Once successfully installed, comprehensive help is available through the Help menu when using the application.

System Requirements

Micro I PC Software requires certain hardware and software components in order to run properly.

PC Requirements

An IBM-compatible PC is required, with hardware that meets or exceeds the following minimum requirements.

Processor: 800 MHz or above

RAM: 256 MB

Free Disk Space: 50 MB

Video: 800x600, 256 colours. It is recommended that a resolution of at least 1280x1024 be used to enable the full benefits of the multi-window interface.

At least one free USB port. (An additional USB port will be required for USB enabled spirometers)

Operating System Requirements

Micro I PC Software will run on the following operating systems:

Windows XP

Windows Vista 32bit

Windows Vista 64bit

Windows Server 2003

Windows Server 2008

Installing Micro I PC Software

Before you begin, please ensure your computer meets the minimum system requirements and the user installing the software has administrative user rights.

Close any other applications that are running.

Insert the installation CD into your CD-Rom drive.

The setup program should launch automatically, displaying the welcome dialog box, click 'Next' to continue. If the setup program does not launch automatically use Windows Explorer to manually select the CD drive then open the file named Micro I PC Software.exe.

The license agreement will be displayed. Please read the complete document and ensure you understand fully before accepting the terms of the license. Click 'Next' to continue.

The select destination location dialog will be displayed showing the directory where the Micro I PC Software will be installed. The default location is:

C:\Program Files\CareFusion\Micro I PC Software

To change this location select 'Browse'.

Click 'Next' to display the Select program manager group dialog. This will be the location where the Micro I PC Software launch icon will be placed under the Start Menu. Either a new location can be entered or an existing location used.

Click 'Next' to display the installation options dialog containing the following option:

Add Micro I PC Software icon to the desktop.


Click 'Next' to start copying the files to your system.

Once the files have finished copying to your system you will be advised that the Micro I PC Software has been successfully installed

and that it would be advisable to restart your PC before using the application. Click 'Finish' to complete the installation process.


Running Micro I PC software

Connect your Micro I to a USB port on the PC using the cable provided. The Micro I will take power from the USB port on the PC. Please be advised that the Micro I requires a High Current USB port, and if connected to a Low Power USB port the device will turn itself off and fail to respond to any key presses until disconnected. Low Power USB ports are commonly found on keyboards and unpowered USB hubs and should not be used. Whilst connected to the PC the Micro I will display:


Do not allow the patient to handle the spirometer during this procedure.

Run the Micro I PC software and the first screen will show:

From this screen you can either change the device settings or adjust the Micro I internal clock. When the change settings options is used the display will change to:


By selecting the 'Custom' option at the top of the screen all of the Micro I features may be tailored to you specific requirements.

Operation

The Micro I is designed to suit a range of applications from the very simplest spirometry test where only a few indices are required to be displayed to more complex operation where deviation from the norm, bronchodilator response, and interpretation of results are required.

The main menu is displayed after the initial configuration and subsequently when the unit is turned on:

Main Menu Overview

Quick Exam

Use this function to take an immediate spirometry measurement with no predicted values or interpretation.

Exam with Predicted

This function requires the entry of the patient's demographics so that predicted values and interpretation may be calculated and displayed.

Post BD Exam

This function allows the post bronchodilator response to be measured. The response is measured with respect to the previously measured baseline obtained using either the Quick Exam or the Exam with Predicted options. The last recorded baseline examination is automatically stored when the unit is turned off and

will be available for a post bronchodilator comparison when the unit is turned on.

Last Results


This option is used to view the results of the last stored examination.

Settings

This option allows the user to adjust various settings including date, time and language and to perform a calibration check.

Quick Exam

After selecting this option the display will change to:

The Micro I may be used directly with a MicroGuard filter or with a disposable cardboard mouthpiece with the adapter supplied. Insert mouthpiece or MicroGuard filter into the mouthpiece holder of the spirometer.

Instruct the patient to inhale as deeply as possible, seal their lips around the mouthpiece and exhale as hard and as fast as possible until no more air can be exhaled.

The spirometry results for that blow are then displayed together with the manoeuvre quality check based upon ATS/ERS guidelines:

14:50		
Current Blow Results		
Base Blow 1		%Pred
FEV1	3.85L	-
FVC	4.55L	-
FEV1/FVC	85%	-
PEF	524L/m	-

Good Blow

Further indices, if configured, and the best results from a sequence of blows may be displayed by pressing the down key.

Each manoeuvre is quality checked for a slow start, abrupt end, short blow, poor effort or cough according to ATS/ERS 2005 guidelines.

Please note that the percentage of predicted value and the interpretation are not available for the Quick Exam option.

Press Enter to display:

11:19

Base testing complete?

Blow Again >

Main Menu >

To repeat the test select Blow Again and press Enter.

Exam with Predicted Values

When this option is selected the following is displayed:

10:56 

Enter Patient Details

Age **25**

Gender Male

Height 175cm

Origin Not-Specified

Use the up and down keys to adjust the patient's age and then press enter. Repeat for gender, height and racial origin. If a mistake is made then simply touch the back key to go back to the previous entry. The racial origin selected applies a percentage drop to the volumetric predicted values to be applied for the patient.

These ethnic corrections for ECCS are as follows:

Population	Percentage drop
Caucasian	0
Non-Specified	0
Polynesian	10
Asian	10
Afro Caribbean	13


If NHANES III predicted values are used then the ethnic origin will be used according to the NHANES III equations.

Press Enter when the required correction has been selected.

Please note that the Micro I may be configured to accept height in inches.

If a mistake is made then simply touch the back key to go back to the previous entry.

When all the patient data has been entered the following screen will be displayed and a spirometry test may be performed:

After performing a test the results are displayed:

```
14:51 ||||  
Current Blow Results  
Base Blow1      %Pred  
FEV1      3.85L  106%  
FVC       4.55L  109%  
FEV1/FVC  85%   101%  
PEF       524L/m 113%  
Predicted:      ECCS  
                Good Blow
```

Each manoeuvre is quality checked for a slow start, abrupt end, short blow, poor effort or cough according to ATS/ERS 2005 guidelines. An asterisk (*) at the end of the line denotes a result below the lower limit of normality.

When a sequence of blows is recorded, the results and the quality check refer to the current blow but the interpretation is based upon the best result of the sequence.

Use the down arrow to see further indices, if configured, and the best results from a sequence of blows:

10:56		
Base Summary		
Best Blow of 2		%Pred
FEV1	4.15L	96%
FVC	4.91L	97%
FEV1/FVC	85%	102%
PEF	566L/m	96%
Predicted:	ECCS	

At any time when the results are displayed another spirometer test may be performed by touching the enter key. If the enter key is touched accidentally then simply press the back key to return to the results screen.

The Micro I may be configured to display the Z score instead of the percentage of predicted values where data from the predicted value sets are available:

14:57		
Current Blow Results		
Base Blow 1		Zscr
FEV1	3.85L	0.4
FVC	4.55L	0.8
FEV1/FVC	85%	0.1
PEF	524L/m	1.1
Predicted:	ECCS	
Good Blow		

The Z score is the number of standard deviations that the results are above (positive Z score) or below (negative Z score) the predicted value.

Post BD Exam

Once satisfactory baseline tests have been recorded a post bronchodilator examination may be performed by selecting the Post BD Exam option from the main menu. When these tests are performed the percentage of predicted value is replaced by the percentage change from the best baseline result for each index:

14:55			
Current Blow Results			
Post Blow 1		%Chg	
FEV1	3.48L	-10%	
FVC	3.88L	-15%	
FEV1/FVC	90%	6%	
PEF	435L/m	-17%	
Predicted:		ECCS	
Possible Cough			

Using the down key will display the other indices and the best results as with the baseline blows.

NLHEP Mode

This mode of operation is obtained by selecting the USA (NLHEP) option from the configuration menu when the unit is first turned on or subsequently from the settings option on the main menu. In this mode the functionality of the Micro I becomes fully compliant with the National Lung Health Education Program requirements.

Only FEV1, FEV6 and the ratio, FEV1/FEV6, are displayed and stored, quality checks are applied to every manoeuvre and a quality grading score for the test session is displayed:

14:54			
Base Summary			
Base Blow2		%Pred	
FEV1	3.8L	103%	
FEV6	4.6L	104%	
FEV1/FEV6	85%	99%	
Predicted:		NHANES III	
NLHEP QC Grade		A	
Good test session			

NLHEP Quality Checks

In order for the interpretation to be displayed a manoeuvre must pass a set of enhanced quality checks. After the patient has performed a poor quality manoeuvre, one of the three following messages will appear:

Message	Criteria	Recommended action
Don't hesitate	Back-extrapolated volume (BEV) greater than 150 ml	The patient should blast out the air more quickly and evenly and without hesitation at the beginning of the manoeuvre
Blast out faster	Time until peak flow (PEFT) greater than 120 msec	The patient must exhale more explosively at the beginning of the manoeuvre
Blow out longer	Expiration time less than 6 seconds or volume accumulation has not dropped below 100 ml per 0.5 seconds	The patient stopped exhaling too early. The patient must exhale until their lungs are completely empty.
Once an acceptable manoeuvre has been performed, the following consistency checks will also be applied to subsequent manoeuvres		
Blast out harder	Peak flow not reproducible. The best previous manoeuvres do not match within 1.0 L/sec indicating that the patient is giving an inconsistent effort.	The patient must give their maximum effort for each manoeuvre
Deeper breath	FEV1 or FEV6 not reproducible. Difference with respect to best test greater than 150 ml	The patient must inhale until their lungs are completely full before each manoeuvre

When two manoeuvres fail either of the consistency checks, then the best individual results of the two are saved (FEV1, FEV6, and PEF individually).

When any of the above messages appear, instruct the patient on how to improve their manoeuvre and demonstrate the correct manoeuvre yourself.

NLHEP QC Grades

The quality of each session is graded according to the following criteria:

QC Grade	Criteria
A	at least two acceptable manoeuvres, with the largest two FEV1's matching better than 100mL and the largest two FEV6's matching better than 100mL
B	at least two acceptable manoeuvres, with the largest two FEV1's matching better than 150mL
C	at least two acceptable manoeuvres, with the largest two FEV1's matching between 200ml and 150mL
D	only one acceptable manoeuvre or two with the largest two FEV1's matching less than 200mL
F	No acceptable manoeuvres

NLHEP Interpretation

The interpretation is performed on the best spirometry results and is based upon the predicted values for the forced expiratory ratio, FEV1/FEV6, and FEV1. If airways obstruction is detected the level of severity is reported in accordance with the NLHEP guidelines.

Switching Off

The unit is switched off by pressing the On/Off button.

The unit can be disconnected from the mains by unplugging the charger from the mains socket or unplugging the USB cable.

Maintenance

Battery Management

The Micro I is powered by a rechargeable battery pack.

The battery voltage is monitored and a warning is displayed on switch on when the battery starts to become exhausted.

Battery Status Icons

Battery nearly exhausted. Recharge as soon as possible to avoid running out of charge.


Battery fully charged.

The battery is not required to hold the internal memory and stored results will not be lost when the battery becomes exhausted.

To recharge the battery connect the Micro I to the charger provided and then plug the charger into a suitable wall socket ensuring that access to the charger is not restricted so that it may easily be removed.

The Micro I will display:

Do not allow the patient to handle the spirometer during this procedure.

The charging icon (🔌) will display to indicate that the device is charging.

Once the charging icon has switched off (up to 5 hours) it will be replaced by the battery fully charged icon. (🔋) Remove the power adapter from the wall socket and from the base of the unit. Micro I is now ready for portable use.


Note: Micro I may be also charged from PC or laptop using the USB cable provided.

Battery Replacement

The lifetime of all rechargeable batteries is limited and the battery pack will need to be replaced after a few years, depending upon usage. As the battery nears the end of its life you will notice that fewer tests can be performed between charges.

For access to the battery the rear of the device requires careful removal.

Clips are moulded into the rear cover in the positions indicated and should be carefully prised undone in the sequence shown.


To aid this process clip 1 has a recess moulded into the device body to allow a small flat bladed screw driver to be carefully slid under the clip to release it.


Once released work the screw driver around to clips 2, 3 and 4.

Finally slide the cover off the bottom of the unit which in turn releases clips 5 and 6

Replace the battery paying particular attention to the orientation of the connector in its socket and of the battery itself in the case. Note the leads coming from the battery should exit from the upper face of the pack when installed.

Replace the rear cover by attaching the securing clips in the reverse order.


Dispose of the waste battery in line with EU Battery Regulations.

Calibration Check

The Micro I is calibrated to read in litres at body temperature, barometric pressure saturated with water vapour (BTPS).

The calibration should remain stable indefinitely, unless the transducer is physically damaged, and the unit should not require re-calibration. However, to ensure the correct functioning of the unit we do recommend that a calibration check is performed periodically.

To perform a calibration check select calibration check option from the setting menu and the display will show:

Connect a 3L syringe to the Micro I with the minimum of adapters and inject the syringe volume into the transducer evenly, without pausing. When complete the display will show:

The acceptable calibration accuracy is $\pm 3.5\%$. If the calibration accuracy is outside of this range then the tick will be replaced by a cross. If this happens then check the syringe and the connections for leaks and repeat the check. If the Micro I is repeatedly outside

the calibration range then the unit will have to be returned to CareFusion for servicing.

Cleaning the Transducer

The transducer requires no routine maintenance or servicing.

To sterilise or clean the transducer it may be removed by means of the following procedure:

- 1 Remove the transducer by gently rotating the transducer anti-clockwise and then pulling from the main body.
- 2 The transducer may now be immersed in warm soapy water for routine cleaning or immersed in cold sterilising solutions e.g. Perasafe for a period not exceeding 10 minutes. (Alcohol and chloride solutions should be avoided.) After cleaning/sterilising, the transducer should be rinsed in distilled water and dried. Perasafe is available from CareFusion in convenient 81g plastic bottles, Cat No. 36-SSC5000A.
- 3 Re-assemble the transducer into the Micro I.

CAUTION: Do not attempt to wash or immerse the spirometer or accessories in water or cleaning fluid, as there are electronic components inside that will be permanently damaged

Servicing

If your unit requires service or repair please see page 212 for contact details.

There are no user serviceable parts in the Micro I.

Trouble Shooting Information

Should you encounter problems operating your spirometer, please consult the table below:

Problem	Possible Cause	Solution
Micro I cannot be switched on	Batteries are flat	Recharge batteries
Every time you switch the instrument on the time is shown as 00:00	The internal battery of Micro I is defective	Contact your dealer.
Micro I is outside of % when conducting the calibration check	There are leaks in the syringe or connections	Check the syringe and connections for leaks

Electromagnetic Compatibility (EMC)

Changes or modifications to the Micro I spirometer that are not expressly approved by Carefusion can cause EMC issues with this or other equipment. The device is designed and tested to comply with applicable regulation regarding EMC and needs to be installed and put into service as per interference exceeding levels specified in EN 50082-1:1992

WARNING: use of portable phones or other radio frequency (RF) emitting equipment near the system may cause unexpected or adverse operation

Safety Designation per IEC 60601-1

Type of protection against electrical shock	Class II
Degree of protection against electrical shock	Type B applied part
Supply connection	100-240 VAC, 50/60Hz
Power Equipment	Adapter and rechargeable internal NiMH battery.
Battery life:	Approximately 30 hours with a fully charged new battery.
Degree of Electrical connection between equipment and Patient	Equipment designed as non-electrical connection to the patient.
Degree of mobility	Transportable
Mode of operation	Continuous

NOTE: When you connect other equipment to the unit, always make sure that the whole combination complies with the international safety standard IEC 60601-1 for medical electrical systems. When connecting to a PC with the supplied USB lead then the PC must be a EN60950-1 compliant SELV limited power source or EN60601-1 compliant.

IMPORTANT: Only use the mains adapter supplied (36-PSU1017 5V DC 1.2A) The adapter contains a transformer. Do not cut off the adapter to replace it with another plug as this causes a hazardous situation.

- The adapter transforms the mains voltage (100-240 Volts) to a safe voltage (5V DC)
- Make sure the adapter does not get wet
- Do not use a damaged adapter
- Always unplug your Micro I before cleaning

WARNING: Do not connect devices that are not specified as part of the system.

WARNING: No modification of this equipment is allowed.

NOTE: If an MPSO (Multiple Portable Socket Outlet) is used with the system, the maximum permitted load should not be exceeded. Do not connect electrical equipment that has not been supplied as part of the system.

Symbols

	Type B device
	In accordance with Directive 93/42/EEC
	Disposal in compliance with WEEE
	Consult the instructions for use
	Caution: consult the accompanying documents
	Date of manufacture
	Manufacturer
	Serial number
	Federal U.S. law restricts this device to sale by or on the order of a physician. (Rx only)

Please Note: Information in this manual is subject to change without notice and does not represent commitment on the part of CareFusion UK 232 Ltd. The software may be used or copied only in accordance with the terms of that agreement. No part of this manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording for any purpose without the written permission of CareFusion UK 232 Ltd.

Specification of the Micro I

Measurements:

Forced Expired Volume in 1 second (FEV1)
Forced Expired Volume in 6 second (FEV6)
Forced Vital Capacity (FVC)
Forced Expiratory Ratio (FEV1/FEV6)
Forced Expiratory Ratio (FEV1/FVC)
Peak Expiratory Flow Rate (PEF),
Mid Expiratory Flow (FEF25-75)
Expiratory Flow at 75% of volume remaining (FEF75)
Expiratory Flow at 25% of volume remaining (FEF25)

Display:	128 x 128 pixel graphic backlit monochrome LCD
Transducer type:	CareFusion Uni-Directional Digital Volume
Accuracy:	To the requirements of the ATS/ERS TASK FORCE: STANDARDISATION OF LUNG FUNCTION TESTING 2005 (Eur Respir J 2005; 26: 319-338 Table 6)
Power supply:	2 X AA size NiMH rechargeable cells
Battery life:	Approximately 30 hours with a fully charged new battery.
Operating current	Less than 90mA
Charging current	Less than 500mA
Dimensions:	162 X 61 X 30mm
Weight:	152g
Operating temperature:	10 to +35° C
Operating humidity:	20% to 80% RH
Storage temperature:	-20 to +70° C
Storage humidity:	30% to 90% RH
Life time	5 years

Accessories

Cat. No.	Description
36-PSA1000	Adult Disposable Mouthpieces (500 per box)
36-SST1250	One-way Safety Mouthpieces (250 per box)
36-SST1000	One-way Safety Mouthpieces (500 per box)
36-MGF1025	Microgard Pulmonary Filter (25 per box)
36-MGF1100	Microgard Pulmonary Filter (100 per box)
36-PSA1200	Paediatric Disposable Mouthpieces (250 per box)
36-PSA1100	Paediatric Adaptor
36-SSC5000A	PeraSafe Sterilising Powder 81g (to make up 5 litres of solution)
36-VOL2104	Nose Clips (pack of 5)
36-SM2125	3 Litre Calibration Syringe
36-MLD1621	Mouthpiece adaptor
36-BAT1043	Battery Pack (NiMH AA 2.4V, 1600mAh)
36-PSU1017	5V DC 1.2A (GlobTek WR9QA1200USBNR-G2267)

PLEASE NOTE: USE ONLY CAREFUSION ACCESSORIES

Micro I Diagnostic Spirometer

Operating Manual – CareFusion USA

Federal (USA) law restricts this device to sale by or on the order of a physician or licensed practitioner.

CareFusion,
22745 Savi Ranch Parkway,
Yorba Linda
CA 92887-4668
USA

II. Introduction – CareFusion USA

The Micro I is a compact, rechargeable battery operated and fully portable diagnostic spirometer. It is accurate to the requirements of the ATS/ERS TASK FORCE: STANDARDISATION OF LUNG FUNCTION TESTING 2005. Its ergonomic and user-friendly design allows diagnostic spirometry measurements, including predicted values and automatic interpretation, to be made quickly and simply.

The Micro I features include:

- Measures FEV₁, FVC(FEV₆), FEV₁/FVC, FEV₁/FEV₆, PEF, FEF₂₅, FEF₇₅, and FEF 25-75
- Fully configurable using software supplied
- ECCS, Asian or NHANES III predicted values
- Variation from norm as percentage or Z-score
- NLHEP, NICE, or the ATS/ERS interpretation
- Spirometry manoeuvre quality checks
- Post bronchodilator comparison
- NLHEP compliant mode

The spirometer uses the CareFusion Digital Volume Transducer, an extremely stable form of volume transducer, which measures expired air directly at B.T.P.S (Body Temperature and Pressure with Saturated water vapour) thus avoiding the inaccuracies of temperature corrections. This transducer is insensitive to the effects of condensation and temperature and avoids the need for individual calibration prior to performing a test.

CareFusion can supply spirometers to fulfil all your diagnostic and monitoring spirometry needs. The Micro II spirometer is the sister spirometer to the Micro I specifically designed for monitoring of asthma and COPD either at home or in the work place.

Package Contents

The Micro I is packaged in a convenient carrying case and comes complete with the following items (Fig.1):

1. Micro I spirometer
2. CareFusion Digital Volume Transducer together with disposable cardboard mouthpieces and instruction manual.
3. Universal power supply (36-PSU1017 5VDC 1.2A)
4. Cardboard mouthpiece adapter
5. USB/Charging cable

Fig 1

Contraindications, Warnings and Cautions

Contraindications: As performing spirometry might be physically demanding, it is contraindicated in patients with recent myocardial infarction. Also extensive exhalation might lead to syncope.

The following terms are used as follows in this manual

Caution: Possibility of injury or serious damage

Warning: conditions or practices that could result in personal injury.

Please Note: Important information for avoiding damage to the instrument or facilitating operation of the instrument.

	CAUTION: Read the manual before use
	WARNING: The instrument is not suitable for use in the presence of explosive or flammable gases, flammable anaesthetic mixtures or in oxygen rich environments.
	CAUTION: Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement.
	CAUTION: Pulmonary filters are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may also increase air resistance and lead to an incorrect measurement.
	CAUTION: Do not allow the patient to handle the spirometer when connected to either the power supply for charging or to a PC when configuring the unit.
	PLEASE NOTE: The product you have purchased should not be disposed of as unsorted waste. Please utilise your local  WEEE collection facilities for the disposal of this product.

Intended Use

The Micro I is used by clinicians and healthcare professionals for the diagnosis of asthma and COPD in adults and children over the age of 3 years.

Environment

Please observe the following precautions:

- Avoid exposing the Micro I to direct sunlight.
- Avoid operating the spirometer in dusty conditions or near to heating appliances or radiators.
- Do not keep the spirometer in a damp place or expose it to extreme temperatures.
- Do not direct the transducer holder towards a strong light source whilst operating the spirometer.

Getting Started

It is recommended that Micro I spirometer is fully charged before use. The power supply is provided with separate UK, USA and European plugs. Connect the required plug to the power supply and plug in to a mains outlet. Connect the Micro I to the power supply using the USB/charging cable and the charging symbol will appear on the screen.

Fully charge for a minimum of 5 hours when used for the first time.

Remove the protective film from the display screen before use.

Configuration

When the unit is first turned on it may be configured for your region. This will set the language, height and weight units, date format and predicted values set appropriate for your region. It will also configure the indices to be displayed and whether percentage predicted or Z-score is to be used to show variation from the norm. However all of these settings may be customised using the PC software supplied.

Turn the unit on by pressing the on/off button located at the top of the device and the following will be displayed:

Use the up and down arrow keys to highlight the required country and then press 'Enter' (↵). This procedure is only required when the unit is first switched on and the selection will be stored for future use.

It is recommended that the PC software be used to adjust these settings, if required. Connect Micro I to the PC using the USB cable supplied. Run the software, and turn on the Micro I. The PC software will detect that the Micro I is connected. Follow the on-screen instructions to configure your unit.

Micro I PC Software

This software allows you to customise the settings on your Micro I spirometer. Once successfully installed, comprehensive help is available through the Help menu when using the application.

System Requirements

Micro I PC Software requires certain hardware and software components in order to run properly.

PC Requirements

An IBM-compatible PC is required, with hardware that meets or exceeds the following minimum requirements.

Processor: 800 MHz or above

RAM: 256 MB

Free Disk Space: 50 MB

Video: 800x600, 256 colours. It is recommended that a resolution of at least 1280x1024 be used to enable the full benefits of the multi-window interface.

At least one free USB port. (An additional USB port will be required for USB enabled spirometers)

Operating System Requirements

Micro I PC Software will run on the following operating systems:

Windows XP

Windows Vista 32bit

Windows Vista 64bit

Windows Server 2003

Windows Server 2008

Installing Micro I PC Software

Before you begin, please ensure your computer meets the minimum system requirements and the user installing the software has administrative user rights.

Close any other applications that are running.

Insert the installation CD into your CD-Rom drive.

The setup program should launch automatically, displaying the welcome dialog box, click 'Next' to continue. If the setup program does not launch automatically use Windows Explorer to manually select the CD drive then open the file named Micro I PC Software.exe.

The license agreement will be displayed. Please read the complete document and ensure you understand fully before accepting the terms of the license. Click 'Next' to continue.

The select destination location dialog will be displayed showing the directory where the Micro I PC Software will be installed. The default location is:

C:\Program Files\CareFusion\Micro I PC Software

To change this location select 'Browse'.

Click 'Next' to display the Select program manager group dialog. This will be the location where the Micro I PC Software launch icon will be placed under the Start Menu. Either a new location can be entered or an existing location used.

Click 'Next' to display the installation options dialog containing the following option:

Add Micro I PC Software icon to the desktop.


Click 'Next' to start copying the files to your system.

Once the files have finished copying to your system you will be advised that the Micro I PC Software has been successfully installed

and that it would be advisable to restart your PC before using the application. Click 'Finish' to complete the installation process.


Running Micro I PC software

Connect your Micro I to a USB port on the PC using the cable provided. The Micro I will take power from the USB port on the PC. Please be advised that the Micro I requires a High Current USB port, and if connected to a Low Power USB port the device will turn itself off and fail to respond to any key presses until disconnected. Low Power USB ports are commonly found on keyboards and unpowered USB hubs and should not be used. Whilst connected to the PC the Micro I will display:


Do not allow the patient to handle the spirometer during this procedure.

Run the Micro I PC software and the first screen will show:

From this screen you can either change the device settings or adjust the Micro I internal clock. When the change settings options is used the display will change to:


By selecting the 'Custom' option at the top of the screen all of the Micro I features may be tailored to you specific requirements.

Operation

The Micro I is designed to suit a range of applications from the very simplest spirometry test where only a few indices are required to be displayed to more complex operation where deviation from the norm, bronchodilator response, and interpretation of results are required.

The main menu is displayed after the initial configuration and subsequently when the unit is turned on:

Main Menu Overview

Quick Exam

Use this function to take an immediate spirometry measurement with no predicted values or interpretation.

Exam with Predicted

This function requires the entry of the patient's demographics so that predicted values and interpretation may be calculated and displayed.

Post BD Exam

This function allows the post bronchodilator response to be measured. The response is measured with respect to the previously measured baseline obtained using either the Quick Exam or the Exam with Predicted options. The last recorded baseline examination is automatically stored when the unit is turned off and

will be available for a post bronchodilator comparison when the unit is turned on.

Last Results


This option is used to view the results of the last stored examination.

Settings

This option allows the user to adjust various settings including date, time and language and to perform a calibration check.

Quick Exam

After selecting this option the display will change to:

The Micro I may be used directly with a MicroGuard filter or with a disposable cardboard mouthpiece with the adapter supplied. Insert mouthpiece or MicroGuard filter into the mouthpiece holder of the spirometer.

Instruct the patient to inhale as deeply as possible, seal their lips around the mouthpiece and exhale as hard and as fast as possible until no more air can be exhaled.

The spirometry results for that blow are then displayed together with the manoeuvre quality check based upon ATS/ERS guidelines:

14:50		
Current Blow Results		
Base Blow 1		%Pred
FEV1	3.85L	-
FVC	4.55L	-
FEV1/FVC	85%	-
PEF	524L/m	-
Good Blow		

Further indices, if configured, and the best results from a sequence of blows may be displayed by pressing the down key.

Each manoeuvre is quality checked for a slow start, abrupt end, short blow, poor effort or cough according to ATS/ERS 2005 guidelines.

Please note that the percentage of predicted value and the interpretation are not available for the Quick Exam option.

Press Enter to display:

11:19		
Base testing complete?		
Blow Again	>	
Main Menu	>	

To repeat the test select Blow Again and press Enter.

Exam with Predicted Values

When this option is selected the following is displayed:

10:56

Enter Patient Details

Age25

GenderMale

Height175cm

OriginNot-Specified

Use the up and down keys to adjust the patient’s age and then press enter. Repeat for gender, height and racial origin. If a mistake is made then simply touch the back key to go back to the previous entry. The racial origin selected applies a percentage drop to the volumetric predicted values to be applied for the patient.

These ethnic corrections for ECCS are as follows:

Population	Percentage drop
Caucasian	0
Non-Specified	0
Polynesian	10
Asian	10
Afro Caribbean	13


If NHANES III predicted values are used then the ethnic origin will be used according to the NHANES III equations.

Press Enter when the required correction has been selected.

Please note that the Micro I may be configured to accept height in inches.

If a mistake is made then simply touch the back key to go back to the previous entry.

When all the patient data has been entered the following screen will be displayed and a spirometry test may be performed:

After performing a test the results are displayed:

```
14:51 ||||  
Current Blow Results  
Base Blow1      %Pred  
FEV1      3.85L  106%  
FVC       4.55L  109%  
FEV1/FVC  85%   101%  
PEF       524L/m 113%  
Predicted:      ECCS  
                Good Blow
```

Each manoeuvre is quality checked for a slow start, abrupt end, short blow, poor effort or cough according to ATS/ERS 2005 guidelines. An asterisk (*) at the end of the line denotes a result below the lower limit of normality.

When a sequence of blows is recorded, the results and the quality check refer to the current blow but the interpretation is based upon the best result of the sequence.

Use the down arrow to see further indices, if configured, and the best results from a sequence of blows:

10:56		
Base Summary		
Best Blow of 2		%Pred
FEV1	4.15L	96%
FVC	4.91L	97%
FEV1/FVC	85%	102%
PEF	566L/m	96%
Predicted:	ECCS	

At any time when the results are displayed another spirometer test may be performed by touching the enter key. If the enter key is touched accidentally then simply press the back key to return to the results screen.

The Micro I may be configured to display the Z score instead of the percentage of predicted values where data from the predicted value sets are available:

14:57		
Current Blow Results		
Base Blow 1		Zscr
FEV1	3.85L	0.4
FVC	4.55L	0.8
FEV1/FVC	85%	0.1
PEF	524L/m	1.1
Predicted:	ECCS	
Good Blow		

The Z score is the number of standard deviations that the results are above (positive Z score) or below (negative Z score) the predicted value.

Post BD Exam

Once satisfactory baseline tests have been recorded a post bronchodilator examination may be performed by selecting the Post BD Exam option from the main menu. When these tests are performed the percentage of predicted value is replaced by the percentage change from the best baseline result for each index:

14:55			
Current Blow Results			
Post Blow 1		%Chg	
FEV1	3.48L	-10%	
FVC	3.88L	-15%	
FEV1/FVC	90%	6%	
PEF	435L/m	-17%	
Predicted:		ECCS	
Possible Cough			

Using the down key will display the other indices and the best results as with the baseline blows.

NLHEP Mode

This mode of operation is obtained by selecting the USA (NLHEP) option from the configuration menu when the unit is first turned on or subsequently from the settings option on the main menu. In this mode the functionality of the Micro I becomes fully compliant with the National Lung Health Education Program requirements.

Only FEV1, FEV6 and the ratio, FEV1/FEV6, are displayed and stored, quality checks are applied to every manoeuvre and a quality grading score for the test session is displayed:

14:54			
Base Summary			
Base Blow2		%Pred	
FEV1	3.8L	103%	
FEV6	4.6L	104%	
FEV1/FEV6 85%		99%	
Predicted:		NHANES III	
NLHEP QC Grade		A	
Good test session			

NLHEP Quality Checks

In order for the interpretation to be displayed a manoeuvre must pass a set of enhanced quality checks. After the patient has performed a poor quality manoeuvre, one of the three following messages will appear:

Message	Criteria	Recommended action
Don't hesitate	Back-extrapolated volume (BEV) greater than 150 ml	The patient should blast out the air more quickly and evenly and without hesitation at the beginning of the manoeuvre
Blast out faster	Time until peak flow (PEFT) greater than 120 msec	The patient must exhale more explosively at the beginning of the manoeuvre
Blow out longer	Expiration time less than 6 seconds or volume accumulation has not dropped below 100 ml per 0.5 seconds	The patient stopped exhaling too early. The patient must exhale until their lungs are completely empty.
Once an acceptable manoeuvre has been performed, the following consistency checks will also be applied to subsequent manoeuvres		
Blast out harder	Peak flow not reproducible. The best previous manoeuvres do not match within 1.0 L/sec indicating that the patient is giving an inconsistent effort.	The patient must give their maximum effort for each manoeuvre
Deeper breath	FEV1 or FEV6 not reproducible. Difference with respect to best test greater than 150 ml	The patient must inhale until their lungs are completely full before each manoeuvre

When two manoeuvres fail either of the consistency checks, then the best individual results of the two are saved (FEV1, FEV6, and PEF individually).

When any of the above messages appear, instruct the patient on how to improve their manoeuvre and demonstrate the correct manoeuvre yourself.

NLHEP QC Grades

The quality of each session is graded according to the following criteria:

QC Grade	Criteria
A	at least two acceptable manoeuvres, with the largest two FEV1's matching better than 100mL and the largest two FEV6's matching better than 100mL
B	at least two acceptable manoeuvres, with the largest two FEV1's matching better than 150mL
C	at least two acceptable manoeuvres, with the largest two FEV1's matching between 200ml and 150mL
D	only one acceptable manoeuvre or two with the largest two FEV1's matching less than 200mL
F	No acceptable manoeuvres

NLHEP Interpretation

The interpretation is performed on the best spirometry results and is based upon the predicted values for the forced expiratory ratio, FEV1/FEV6, and FEV1. If airways obstruction is detected the level of severity is reported in accordance with the NLHEP guidelines.

Switching Off

The unit is switched off by pressing the On/Off button.

The unit can be disconnected from the mains by unplugging the charger from the mains socket or unplugging the USB cable.

Maintenance

Battery Management

The Micro I is powered by a rechargeable battery pack.

The battery voltage is monitored and a warning is displayed on switch on when the battery starts to become exhausted.

Battery Status Icons

Battery nearly exhausted. Recharge as soon as possible to avoid running out of charge.


Battery fully charged.

The battery is not required to hold the internal memory and stored results will not be lost when the battery becomes exhausted.

To recharge the battery connect the Micro I to the charger provided and then plug the charger into a suitable wall socket ensuring that access to the charger is not restricted so that it may easily be removed.

The Micro I will display:

Do not allow the patient to handle the spirometer during this procedure.

The charging icon (🔌) will display to indicate that the device is charging.

Once the charging icon has switched off (up to 5 hours) it will be replaced by the battery fully charged icon. (🔋) Remove the power adapter from the wall socket and from the base of the unit. Micro I is now ready for portable use.


Note: Micro I may be also charged from PC or laptop using the USB cable provided.

Battery Replacement

The lifetime of all rechargeable batteries is limited and the battery pack will need to be replaced after a few years, depending upon usage. As the battery nears the end of its life you will notice that fewer tests can be performed between charges.

For access to the battery the rear of the device requires careful removal.

Clips are moulded into the rear cover in the positions indicated and should be carefully prised undone in the sequence shown.


To aid this process clip 1 has a recess moulded into the device body to allow a small flat bladed screw driver to be carefully slid under the clip to release it.


Once released work the screw driver around to clips 2, 3 and 4.

Finally slide the cover off the bottom of the unit which in turn releases clips 5 and 6

Replace the battery paying particular attention to the orientation of the connector in its socket and of the battery itself in the case. Note the leads coming from the battery should exit from the upper face of the pack when installed.

Replace the rear cover by attaching the securing clips in the reverse order.


Dispose of the waste battery in line with EU Battery Regulations.

Calibration Check

The Micro I is calibrated to read in litres at body temperature, barometric pressure saturated with water vapour (BTPS).

The calibration should remain stable indefinitely, unless the transducer is physically damaged, and the unit should not require re-calibration. However, to ensure the correct functioning of the unit we do recommend that a calibration check is performed periodically.

To perform a calibration check select calibration check option from the setting menu and the display will show:

Connect a 3L syringe to the Micro I with the minimum of adapters and inject the syringe volume into the transducer evenly, without pausing. When complete the display will show:

The acceptable calibration accuracy is $\pm 3.5\%$. If the calibration accuracy is outside of this range then the tick will be replaced by a cross. If this happens then check the syringe and the connections for leaks and repeat the check. If the Micro I is repeatedly outside

the calibration range then the unit will have to be returned to CareFusion for servicing.

Cleaning the Transducer

The transducer requires no routine maintenance or servicing.

To sterilise or clean the transducer it may be removed by means of the following procedure:

1. Remove the transducer by gently rotating the transducer anti-clockwise and then pulling from the main body.
2. The transducer may now be immersed in warm soapy water for routine cleaning or immersed in cold sterilising solutions e.g. Perasafe for a period not exceeding 10 minutes. (Alcohol and chloride solutions should be avoided.) After cleaning/sterilising, the transducer should be rinsed in distilled water and dried. Perasafe is available from CareFusion in convenient 81g plastic bottles, Cat No. 36-SSC5000A.
3. Re-assemble the transducer into the Micro I.

CAUTION: Do not attempt to wash or immerse the spirometer or accessories in water or cleaning fluid, as there are electronic components inside that will be permanently damaged

Servicing

If your unit requires service or repair please see page 212 for contact details.

There are no user serviceable parts in the Micro I.

Trouble Shooting Information

Should you encounter problems operating your spirometer, please consult the table below:

Problem	Possible Cause	Solution
Micro I cannot be switched on	Batteries are flat	Recharge batteries
Every time you switch the instrument on the time is shown as 00:00	The internal battery of Micro I is defective	Contact your dealer.
Micro I is outside of % when conducting the calibration check	There are leaks in the syringe or connections	Check the syringe and connections for leaks

Electromagnetic Compatibility (EMC)

Changes or modifications to the Micro I spirometer that are not expressly approved by Carefusion can cause EMC issues with this or other equipment. The device is designed and tested to comply with applicable regulation regarding EMC and needs to be installed and put into service as per interference exceeding levels specified in EN 50082-1:1992

WARNING: use of portable phones or other radio frequency (RF) emitting equipment near the system may cause unexpected or adverse operation

Safety Designation per IEC 60601-1

Type of protection against electrical shock	Class II
Degree of protection against electrical shock	Type B applied part
Supply connection	100-240 VAC, 50/60Hz
Power Equipment	Adapter and rechargeable internal NiMH battery.
Battery life:	Approximately 30 hours with a fully charged new battery.
Degree of Electrical connection between equipment and Patient	Equipment designed as non-electrical connection to the patient.
Degree of mobility	Transportable
Mode of operation	Continuous

NOTE: When you connect other equipment to the unit, always make sure that the whole combination complies with the international safety standard IEC 60601-1 for medical electrical systems. When connecting to a PC with the supplied USB lead then the PC must be a EN60950-1 compliant SELV limited power source or EN60601-1 compliant.

IMPORTANT: Only use the mains adapter supplied (36-PSU1017 5V DC 1.2A) The adapter contains a transformer. Do not cut off the adapter to replace it with another plug as this causes a hazardous situation.

- The adapter transforms the mains voltage (100-240 Volts) to a safe voltage (5V DC)
- Make sure the adapter does not get wet
- Do not use a damaged adapter
- Always unplug your Micro I before cleaning

WARNING: Do not connect devices that are not specified as part of the system.

WARNING: No modification of this equipment is allowed.

NOTE: If an MPSO (Multiple Portable Socket Outlet) is used with the system, the maximum permitted load should not be exceeded. Do not connect electrical equipment that has not been supplied as part of the system.

Symbols

	Type B device
	In accordance with Directive 93/42/EEC
	Disposal in compliance with WEEE
	Consult the instructions for use
	Caution: consult the accompanying documents
	Date of manufacture
	Manufacturer
	Serial number
	Federal U.S. law restricts this device to sale by or on the order of a physician. (Rx only)

Please Note: Information in this manual is subject to change without notice and does not represent commitment on the part of CareFusion UK 232 Ltd. The software may be used or copied only in accordance with the terms of that agreement. No part of this manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording for any purpose without the written permission of CareFusion UK 232 Ltd.

Specification of the Micro I

Measurements:

Forced Expired Volume in 1 second (FEV1)
Forced Expired Volume in 6 second (FEV6)
Forced Vital Capacity (FVC)
Forced Expiratory Ratio (FEV1/FEV6)
Forced Expiratory Ratio (FEV1/FVC)
Peak Expiratory Flow Rate (PEF),
Mid Expiratory Flow (FEF25-75)
Expiratory Flow at 75% of volume remaining (FEF75)
Expiratory Flow at 25% of volume remaining (FEF25)

Display:	128 x 128 pixel graphic backlit monochrome LCD
Transducer type:	CareFusion Uni-Directional Digital Volume
Accuracy:	To the requirements of the ATS/ERS TASK FORCE: STANDARDISATION OF LUNG FUNCTION TESTING 2005 (Eur Respir J 2005; 26: 319-338 Table 6)
Power supply:	2 X AA size NiMH rechargeable cells
Battery life:	Approximately 30 hours with a fully charged new battery.
Operating current	Less than 90mA
Charging current	Less than 500mA
Dimensions:	162 X 61 X 30mm
Weight:	152g
Operating temperature:	10 to +35° C
Operating humidity:	20% to 80% RH
Storage temperature:	-20 to +70° C
Storage humidity:	30% to 90% RH
Life time	5 years

Accessories

Cat. No.	Description
36-PSA1000	Adult Disposable Mouthpieces (500 per box)
36-SST1250	One-way Safety Mouthpieces (250 per box)
36-SST1000	One-way Safety Mouthpieces (500 per box)
36-MGF1025	Microgard Pulmonary Filter (25 per box)
36-MGF1100	Microgard Pulmonary Filter (100 per box)
36-PSA1200	Paediatric Disposable Mouthpieces (250 per box)
36-PSA1100	Paediatric Adaptor
36-SSC5000A	PeraSafe Sterilising Powder 81g (to make up 5 litres of solution)
36-VOL2104	Nose Clips (pack of 5)
36-SM2125	3 Litre Calibration Syringe
36-MLD1621	Mouthpiece adaptor
36-BAT1043	Battery Pack (NiMH AA 2.4V, 1600mAh)
36-PSU1017	5V DC 1.2A (GlobTek WR9QA1200USBNR-G2267)

For further information or to place an order for Disposables/ Supporting products please contact CAREFUSION, your local distributor or view our website www.carefusion.com/micromedical

For US Customer Care: Toll Free 1-800-231-2466, phone 1-714-283-2228

PLEASE NOTE: USE ONLY CAREFUSION ACCESSORIES

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
The screenshot shows the main menu of the Micro I device. At the top left, the time '08:52' is displayed. At the top right, there is a battery status icon. The menu items are listed below, with 'Schnell-Untersuchung' highlighted in a black box. The items are: 'Schnell-Untersuchung', 'Untersuch. mit Sollwerte :', 'Untersuchung nach BD :', 'Letzte Ergebnisse :', and 'Einstellungen :'. Each item is followed by a colon and a space.

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Type transducer: CareFusion Uni-Directional Digital Volume

Nauwkeurigheid: Voldoet aan de eisen van de ATS/ERS TASK FORCE: STANDARDISATION OF LUNG FUNCTION TESTING 2005 (Eur Respir J 2005; 26: 319-338 tabel 6)

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Customer contact information

UK Customers only

For all Sales Order processing for products, training and spare parts,
Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd
UK Customer Service & Support
The Crescent
Jays Close
Basingstoke
RG22 4BS

Customer Service Sales Enquiries:

Telephone: Freephone 0800 151 3584

Telephone: 01256 388550

Email: micro.uksales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: Freephone 0800 151 3585

Email: micro.ukservice@carefusion.com

Technical Support Enquiries:

Telephone: 01256 388551

Email: support.rt.eu@carefusion.com

International customers only

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

Carefusion Germany 234 GmbH
Customer Service & Support International
Leibnizstrasse 7
97204 Hoechberg
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670

Email: micro.internationalsales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867

Email: support.admin.eu@carefusion.com

Technical Support Enquiries:

Telephone: 0049 931 4972 127

Email: support.rt.eu@carefusion.com

CareFusion UK 232 Ltd.
Quayside
Chatham Maritime
Kent ME4 4QY
U.K.

CE
0086

JDE no. 36-MAN1290
Drg no. 122-14WW (PART 1)
Issue 1.0
March 2011
© CareFusion 2011